

HEART INFUSION (PORCINE)

CAT. N°: 303

DESCRIPTION:

Heart Infusion (Porcine) is a specially formulated preparation of dried infusion from pork hearts supplied as a light, fine powder. It is a suitable substitute for heart infusion (bovine) in select microbiological culture media.

CHEMICAL CHARACTERISTICS

SPECIFICATIONS

TYPICAL ANALYSIS

CHEMICAL CHARACTERISTICS	SPECIFICATIONS	TYPICAL ANALYSIS
Total Nitrogen (TN)	Minimum 10.0%	11.7%
Loss on drying	Maximum 6.0%	4.0%
Ash	Maximum 15.0%	10.50%
pH (2% solution)	6.5 - 7.5	7.1

ELEMENTAL PROFILE

Calcium	0.018%
Magnesium	0.012%
Potassium	2.36%
Sodium	3.96%

GROWTH SUPPORTING PROPERTIES

Peptone agar	Satisfactory
--------------	--------------

MICROBIOLOGICAL ANALYSIS

Standard plate count	Less than 5000 CFU/g
Yeasts and molds	Less than 100 CFU/g
Coliforms	Negative
Salmonella	Negative

AMINO ACIDS

TOTAL g/100 g

Alanine	3.32
Arginine	3.08
Aspartic acid	8.77
Cystine	0.54
Glutamic acid	14.77
Glycine	2.55
Histidine	0.05
Isoleucine	3.39
Leucine	5.83
Lysine	5.43
Methionine	1.65
Phenylalanine	3.15
Proline	5.84
Serine	3.81
Threonine	3.33
Tryptophan	0.91
Tyrosine	1.44
Valine	4.24