

GRANULAR HEMOGLOBIN

CAT. N°: 012

DESCRIPTION:

This product is obtained from desiccated bovine blood. It is used in preparation of bacteriological culture media as sources of nutrients for cultivation and isolation of fastidious microorganisms, especially *Neisseria* spp and *Haemophilus* spp.

CHEMICAL CHARACTERISTICS

SPECIFICATIONS

TYPICAL ANALYSIS

CHEMICAL CHARACTERISTICS	SPECIFICATIONS	TYPICAL ANALYSIS
Amino Nitrogen (AN)	Maximum 5.0%	1.0%
Total Nitrogen (TN)	7.8 - 8.6	8.2

GROWTH SUPPORTING PROPERTIES

Chocolate agar	Satisfactory
----------------	--------------